

# BELIRINSE 444

BELIRINSE 444 is a mildly acidic rinse agent used in automated wash equipment to promote accelerated spot free drying of utensils, glassware, accessories and carts.

Contains citric acid, chelating agent & surfactants and is compatible with stainless steel, anodized aluminum, zinc, non-ferrous metals, and other delicate materials

## FEATURES

Contains rinsing agent  
Contains surfactants  
Highly concentrated  
Mild  
Low foaming  
Biodegradable

## BENEFITS

Rinses spot free in all water hardness conditions  
Great wetting ability for free sheeting action  
Low cost per use  
Effectively rinses instruments without harm  
Effective with all automated washers  
Non-polluting, ecology minded

## DIRECTIONS FOR USE

Use 3.5 ml to 5 ml per liter of water or 1/2 ounce to 2/3 ounce per gallon of water.  
Dilutions may vary due to water conditions and type of equipment.  
Contact your Belimed Sales Consultant for recommended usages.

## PHYSICAL PROPERTIES

Appearance:	Clear yellow liquid
Fragrance:	None
pH 1% @ 20°C:	3.19
SG @ 20°C	1.01
Foam:	Low
Available Alkalinity:	None

Available in 5, 15, 30, and 55 gallon containers