

# BELICHEM 900

BELICHEM 900 is a neutral liquid detergent for use in automated wash equipment for the removal of moderate soil contamination on carts and accessories.

Contains surfactant, builders and chelating agent for good wetting ability and detergent penetration. Safe to use on stainless steel, anodized aluminum, zinc, rubber, latex, and other delicate metals.

## FEATURES

Contains chelating agent  
Contains surfactants  
Highly concentrated  
Mild  
Low foaming  
Biodegradable

## BENEFITS

Hard water tolerant  
Great wetting ability for optimal penetration  
Low cost per use  
Effectively rinses equipment without harm  
Effective with all automated washers  
Non-polluting, ecology minded

## DIRECTIONS FOR USE

Use 3.5 ml to 5 ml per liter of water or 1/2 ounce to 2/3 ounce per gallon of water.  
Dilutions may vary due to water conditions and type of equipment.  
Contact your Belimed Sales Consultant for recommended usages.

## PHYSICAL PROPERTIES

Appearance:	Clear blue liquid
Fragrance:	Fresh, floral
pH 1% @ 20°C:	8.8
SG @ 20°C	1.06
Foam:	Low
Available Alkalinity:	Negligible

Available in 5, 15, 30, and 55 gallon containers